

VIA Anti-Gravity Treadmill™

The AlterG® Via Anti-Gravity Treadmill changes what's possible in rehabilitation and training. With unmatched precision and control, the Via combines revolutionary Differential Air Pressure (DAP) technology, an intuitive touchscreen interface, real-time gait analytics, and live video monitoring – all to enable the ultimate in pain-free and fall-safe body weight supported walking, running, and closed-kinetic chain exercise. Now you can help diverse ortho and neuro patients, active seniors, and athletes to gain mobility, recover smarter, and achieve performance goals sooner.

NEW DESIGN

Accessible Design

Snap-up, magnetized air chamber closure lays flat for ease of entry, whilst a new window design provides critical visual accessibility in a robust, functional design.

Smooth Setup

Open-back cockpit design with light and easy height adjustment leverages inseam set-point measurements and smooth-set locking mechanism to make entry and setup simple.

Easy In, Easy Out

New, low step-up height ensures the easiest possible in-and-out for users.

Enhanced Data

AlterG Stride Smart provides valuable insight into gait and session data, allowing user-progress evaluation with essential data metrics.



Functional Safety

AlterG's patented DAP technology reduces gravitational forces and provides a fall-safe environment, whilst supporting unimpeded gait mechanics and balance.



VIA X
Anti-Gravity Treadmill

ANTI-GRAVITY PLATFORM



SPECIFICATIONS

Body Weight Adjustment	100% down to as little as 20% in 1% increments		
Display	15"	15"	21"
Treadmill	Flat belt	Flat belt	Slat belt
Speed	0 mph to 12 mph	-3 mph to 12 mph	-10 mph to 18 mph
Starting Speed	0.1 mph	0.1 mph	0.1 mph
X Package (High Speed)*	Optional (0 mph -15 mph, 0.5 mph Starting Speed, Fit X Graphic)	Optional (-5 mph to 15 mph, 0.5 mph Starting Speed, Via X Graphic)	-
Reverse (-3 mph)	-	Included	Included
Incline (0% - 15%)	-	Included	Included
Live Video Monitoring*	-	Included	Included
Stride Smart	-	Included	Included
AlterG Assistant	-	Included	Included
User Weight	85 lbs to 400 lbs	85 lbs to 400 lbs	85 lbs to 400 lbs

*Service call required for field upgrade(s).

Learn more at alterg.com or call +1.510.270.5900.