

h/p/cosmos® Laufband-Ergom. Biomechanik	Artikelnummer	Lauffläche L/B	Geschwindigkeit	Steigung	Versorgung	Display	CE
pluto® lt	cos30027-01va01	150 / 50 cm	0 ... 18 km/h	0 ... 20 %	200 ... 240 Volt AC 1~ / 16 Amp.	nein	CE
pluto®	cos30026-01va01	150 / 50 cm	0 ... 18 km/h	0 ... 20 %	200 ... 240 Volt AC 1~ / 16 Amp.	ja	CE
pluto® ef	cos30027-01va13	150 / 50 cm	0 ... 18 km/h	0 ... 20 %	200 ... 240 Volt AC 1~ / 16 Amp.	ja	CE
stratos®	cos30000-02va09	150 / 50 cm	0 ... 22 km/h	0 %	200 ... 240 Volt AC 1~ / 16 Amp.	ja	CE
mercury®	cos30000-02va01	150 / 50 cm	0 ... 22 km/h	0 ... 25 %	200 ... 240 Volt AC 1~ / 16 Amp.	ja	CE
quasar®, mit MCU5 (6 Displays)	cos30003va16	170 / 65 cm	0 ... 25 km/h	0 ... 28 %	230 Volt AC 1~ / 16 Amp.	ja	CE
stellar®, mit MCU6 (Touch-Display)	cos30003-01va05	170 / 65 cm	0 ... 25 km/h	0 %	230 Volt AC 1~ / 16 Amp.	Touch	CE
quasar® lt, mit MCU6	cos30003-01va03	170 / 65 cm	0 ... 25 km/h	0 ... 28 %	230 Volt AC 1~ / 16 Amp.	nein	CE
quasar®, mit MCU6 (Touch-Display)	cos30003-01va01	170 / 65 cm	0 ... 25 km/h	0 ... 28 %	230 Volt AC 1~ / 16 Amp.	Touch	CE
stellar® sport 190/65, mit MCU6 (Touch)	cos30004-01va05	190 / 65 cm	0 ... 40 km/h	0 %	230 Volt AC 1~ / 16 Amp.	Touch	CE
pulsar® lt sport, mit MCU6	cos30004-01va03	190 / 65 cm	0 ... 40 km/h	-25 ... +25 %	230 Volt AC 1~ / 16 Amp.	nein	CE
pulsar® sport, mit MCU6 (Touch)	cos30004-01va01	190 / 65 cm	0 ... 40 km/h	-25 ... +25 %	230 Volt AC 1~ / 16 Amp.	Touch	CE
saturn® 450/300 rs sport, mit MCU6 (Touch)	cos30013-01va02	450 / 300 cm	0 ... 40 km/h	-4 ... +25 %	400 Volt AC 3~ / 64 Amp.	Touch	CE

h/p/cosmos® Laufband-Ergom. Biomechanik	Artikelnummer	Lauffläche L/B	Geschwindigkeit	Steigung	Versorgung	Display	CE
gaitway® 1d 3G 150/50 Biomechanik-, inkl. 1-komponenten (Fz) vertikale Bodenreaktionskraft-Messung + Tandem Kraftmessplatten	cos30002-01	150 / 50 cm	0 ... 18 km/h	0 ... 20 %	200 ... 240 Volt AC 1~ / 16 Amp.	ja	CE
stratos®	cos30000-02va09	150 / 50 cm	0 ... 22 km/h	optional	200 ... 240 Volt AC 1~ / 16 Amp.	ja	CE
+ gaitway 3d Biomechanik-Upgrade + 3 Komponenten (Fx,y,z) Kraftmessung	cos102999_150-50_G6	150 / 50 cm	0 ... 22 km/h	optional	200 ... 240 Volt AC 1~ / 16 Amp.	ja	CE
stellar®, mit MCU6 (Touch-Display)	cos30003-01va05	170 / 65 cm	0 ... 25 km/h	optional	230 Volt AC 1~ / 16 Amp.	Touch	CE
+ gaitway 3d Biomechanik-Upgrade + 3 Komponenten (Fx,y,z) Kraftmessung	cos102999_170-65_MCU6	170 / 65 cm	0 ... 25 km/h	optional	230 Volt AC 1~ / 16 Amp.	Touch	CE
stellar® sport 190/65, mit MCU6 (Touch-Display)	cos30004-01va05	190 / 65 cm	0 ... 40 km/h	optional	230 Volt AC 1~ / 16 Amp.	Touch	CE
+ gaitway 3d Biomechanik-Upgrade + 3 Komponenten (Fx,y,z) Kraftmessung	cos102999_190-65	190 / 65 cm	0 ... 40 km/h	optional	230 Volt AC 1~ / 16 Amp.	Touch	CE
pulsar® 3p, mit MCU5 (6 Displays)	cos30004va04	190 / 65 cm	0 ... 40 km/h	optional	400 Volt AC 3~ / 16 Amp.	ja	CE
+ gaitway 3d Biomechanik-Upgrade + 3 Komponenten (Fx,y,z) Kraftmessung	cos102999_190-65	190 / 65 cm	0 ... 40 km/h	optional	400 Volt AC 3~ / 16 Amp.	ja	CE