

WIITEST

This platform is designed to assess muscle strength levels in real time. Its modularity allows it to perform different tests that are scientifically recognised and integrated into injury prevention and performance monitoring protocols.

Dimensions :

105 x 110 x 119 cm
(Depth x Width x Height)

Weight :

Structure 68kg / Platform 22kg

Other features :

Removable platform /
Easy to move
(wheels and handles)



- **Rotation adjustment :**
13 positions / Every 15°
- **Horizontal adjustment :**
11 positions / Every 3cm (per side)
- **Vertical adjustment :**
22 positions / Every 3.5cm



KINVENT
MEASURE. MOVE. PROGRESS.

